

INNER EAR MRI PROTOCOL

Brain MRI

- Axial FLAIR
- Axial T1 SE
- Axial T2 TSE
- Coronal T2 TSE
- Axial DWI
- Axial T2 gradient if necessary

IAM:

- Axial T2TSE 3D – thin slices
- Coronal T2 TSE 3D – thin slices
- Sagittal T2 TSE 3D – (left & right) – planning perpendicular on IAM – also thin slices
- Axial T1 TFE 3D – thin slices

AFTER CONTRAST

IAM

- Axial T1 TFE 3D – thin slices

Brain

- Axial, coronal and sagittal T1 SE
- This procedure takes about 50 minutes